HGX-200 Cold Vapor & Hydride Generator Safety Manual

Original Instructions

Document part number 480209 rev0, Printed in USA

About This Manual

This manual contains instructions to safely operate the HGX-200. This document helps the product comply with the Machinery Directive (2006/42/EC, Europe) and the Equipment and Product Safety Act (ProdSG, Germany).

See the English language *HGX-200 Cold Vapor and Hydride Generator Operator's Manual* for additional instructions on setting up and using the product.

Operator Qualifications

Only analytical chemists and lab technicians should use the HGX-200. Operators must understand laboratory safety, chemical sample handling, and all equipment used with the HGX-200.

WARNING

CHEMICAL INJURY HAZARD

Operators must be trained in laboratory safety. Operators must know the hazards of the chemicals in the laboratory, and take the appropriate precautions. Laboratory chemicals may cause serious injury.

Review this product and related documentation to familiarize with safety markings and instructions before you operate the HGX-200.

The HGX-200 does not contain mechanical or electrical components. The primary hazards associated with the product are chemical injury hazards.

This product is manufactured by Teledyne CETAC Technologies, 14306 Industrial Road, Omaha, Nebraska 68144, USA.

Safety Notices

WARNING

INJURY HAZARD FROM IMPROPER USE

If the equipment is used in a manner not specified by Teledyne CETAC Technologies, the protection provided with the equipment may be impaired.

Repair or service that this not covered in the manual should only be performed by qualified personnel.

WARNING**CHEMICAL HAZARDS**

Learn about the chemicals which will be used in and near the equipment, and observe the necessary precautions. Always wear appropriate personal protective equipment, including face shield or goggles, gloves, and a laboratory coat. Ensure that engineering controls (such as ventilation systems) are operating properly. Do not mix waste borohydride with acids.

WARNING**INJURY HAZARD FROM HYDRIDE EMISSIONS**

Please refer to the MSDS data sheets for metal hydrides which may be present in the laboratory, such as AsH_3 , SbH_3 , H_2Se , $\text{Hg}(0)$, or H_2Te . Provide adequate ventilation to prevent exposure to hazardous vapors or gases.

Note regarding Poisonous Hydride Emissions: Since hydride generation is typically used to analyze trace levels (ppb levels) of potentially poisonous metalloids and metals, it is highly unlikely that any dangerous accumulation of such gases will occur in a lab environment at the typical sample consumption rates of only a few mL/min or ng/min. The emissions of the plasma gases do not contain any hydrides anymore and are no different from normal use of an ICP. Higher levels of metal hydrides may occur if waste borohydride reagent is accidentally mixed with waste solutions.

WARNING**INJURY HAZARD FROM GLASSWARE**

Handling of glassware: Like with any glassware, it is important to avoid stressing joints or applying pressure on connectors. Do not use tools to tighten the fittings.

CAUTION**AVOID DRAIN OVERFLOW:**

This is the most serious risk in terms of creating problems for the ICP/ICPMS equipment. Because the sample and reagent flow rates are quite high (a few mL/min) and the internal volume of the HGX-200 is limited, it takes only a few minutes of improper drainage to cause the unit to overflow. This can extinguish the ICP plasma and could lead to aggressive solutions flowing into the ICP torch box area. This is particularly bad for units with the torch box mounted underneath the RF load coil like on the THERMO ELEMENT2.

Every time you use the HGX-200, make absolutely sure that the drain is pumping sufficiently. Change the pump tubing often (PHARMED® at least every 40 hours of operation, PVC/Tygon® every 20 hours). This is especially true for users that use this equipment only occasionally.

Equipment Startup Procedure

FIRST THING EVERY TIME:

Before you even think about starting the ICP, run the peristaltic pump with sample and reagents flowing, and check the drainage!

Returning the Product for Service

Contact the manufacturer for instructions before returning the product for service.

Do not return products which are contaminated by radioactive materials, infectious agents, or other materials constituting health hazards to CETAC employees.

Explanation of Caution and Warning Notices

Warning symbol marked on equipment. This symbol means “Attention! Refer to the manual.”

WARNING

The **WARNING** notice calls attention to conditions that could result in personal injury. Do not proceed beyond a **WARNING** notice until the indicated conditions are fully understood.

CAUTION

The **CAUTION** notice calls attention to conditions that could result in equipment damage, loss of data, or inaccurate data.

Explanation of Regulatory Marks

The CE mark is a registered trademark of the European Community. This CE mark, if present on the product, shows that the product complies with all the relevant European Legal Directives.

PharMed and Tygon are registered trademarks of Saint-Gobain Performance Plastics.

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