

ADX-500 ONLINE AUTODILUTION ACCESSORY



ACCESSORIES AND PARTS CATALOG

Accessories and Supplies



Z-axis drive assembly (pull-pull style, dual probe) for units with serial numbers >080101ASX and <040248ASX.

SP6015

Sample probe not included.



Z-axis drive assembly (dual probe) for units with serial numbers < 080101ASX and >040249ASX.

SP5956

Sample probe not included.



Sample probe, 0.9 mm ID x 70" made with Ultem brand polyetherimide resin (yellow band)

SP5950



Sample probe, 0.9mm ID x 70", carbon fiber reinforced (yellow band)

SP5950C

All wetted surfaces (liquid contact) of both the carbon fiber and Ultem probes are inert and identical (PTFE). Color bands indicate inside diameters and are the same for both probe types.

ADX-500 ACCESSORIES AND PARTS



3-way valve
For serial numbers < 080601ADX
No syringe is available to go with this valve

SP5951



3-way valve
For serial numbers \geq 080601ADX to < 081101ADX
Goes with syringe SP5952)

SP6529



3-way valve
For serial numbers \geq 081101ADX

SP6885




5 mL syringe
For serial numbers \geq 080601ADX to
< 081101ADX

SP5952



5 mL syringe
For serial numbers \geq 081101ADX

SP6886



Z-axis stepper motor assembly

SP5965

ADX-500 ACCESSORIES AND PARTS




Z-axis rotor disc assembly

SP6596



Z-Axis Rotor Clamp and Screw Kit

SP5964



Switch / power input kit (single pin input)

SP5941

For serial numbers < 090801ADX



Switch / Power input kit (6 pin input)

SP6104

For serial numbers \geq 090801ADX



I/O board

SP5963

For serial numbers < 081101ADX



Main processor board (includes firmware)

SP5962

For serial numbers < 081101ADX

ADX-500 ACCESSORIES AND PARTS

 A green printed circuit board (PCB) with various electronic components, including a microcontroller, capacitors, and connectors. It is labeled as an ADX Control board with USB.	ADX Control board w/USB (includes firmware) For serial numbers \geq 081101ADX	SP6888
 A black, rectangular syringe pump with a white syringe attached to the top. It has a circular dial on the front and a 3-way valve on top.	Syringe pump (includes 3-way valve) For serial numbers $<$ 081101ADX	SP5961
 A black, rectangular syringe pump with a white syringe attached to the top. It has a circular dial on the front and a 3-way valve on top.	Syringe pump (includes 3-way valve) For serial numbers \geq 081101ADX	SP6887
 A black, rectangular power supply unit with a power cord and a connector cable.	24-Volt Power Supply For serial numbers \geq 090801ADX	SP6103
 A black, rectangular power supply unit with a power cord and a connector cable.	24 volt power supply For serial numbers $<$ 090801ADX	SP5526

ADX-500 ACCESSORIES AND PARTS



Serial cable (Female/Female)

SP5341



USB Cable (3 Meters)

SP6375



Diluent bottle (1 L)

SP5957



Diluent uptake tubing (Tygon®)

SP6011

ADX-500 ACCESSORIES AND PARTS

Ordering

To order, contact:

Phone (Toll-Free, USA): 800-369-2822
Phone (Outside USA): +1-402-733-2829
Fax: +1-402-733-5292
e-mail: cetacsales@teledyne.com
www.cetac.com

Teledyne CETAC Technologies
14306 Industrial Road
Omaha, Nebraska 68144, USA

For questions about service and repair, call us at the numbers above or e-mail
cetacservic@teledyne.com

Trademark Acknowledgements: Tygon is a registered trademark of Saint-Gobain Performance Plastics.