

SDS-550 OFFLINE SAMPLE DILUTION SYSTEM ACCESSORIES AND SUPPLIES CATALOG



Teledyne CETAC Technologies
14306 Industrial Road
Omaha, Nebraska 68144, USA
Phone: (402) 733-2829
Fax: (402) 733-5292
e-mail: cetacsales@teledyne.com
www.teledynecetac.com

August 2018

PharMed and Tygon are registered trademarks of
Saint Gobain Performance Plastics


SDS-550 SPARE PARTS

IMAGE	DESCRIPTION	PART #
	<p>Sample Loop 10mL</p>	<p>SP5487</p>
	<p>Ferrule / Fittings Kit</p>	<p>SP6434</p>
	<p>Sample Probe, 1.0mm ID (2 blue bands) {Thermo UV}</p> <p>Sample Probe, 1.0mm ID (2 blue bands) {Mercury}</p> <p>Carbon Probe, 1.0mm ID (2 blue bands) {Thermo UV}</p> <p>Carbon Probe 1.0mm ID (2 blue bands) {Mercury}</p>	<p>SP6013</p> <p>SP6033</p> <p>SP6013C</p> <p>SP6033C</p>
	<p>Solenoid Control Cable</p>	<p>SP6418</p>

SDS-550 SPARE PARTS

IMAGE	DESCRIPTION	PART #
	24-Volt Power Supply	SP5526
	Switch / Power Input	SP6104
	XLP-3000 Syringe Pump	SP6372
	<p>XLP-3000 10mL Syringe Units \geq 020701SDS</p> <p>XL-3000 10mL Syringe Units $<$ 020701SDS</p>	<p>SP6114A</p> <p>SP6114</p>

SDS-550 SPARE PARTS

IMAGE	DESCRIPTION	PART #
	<p>3-Way Valve Assy. For the XLP-3000 Syringe Units \geq 020701SDS (or have upgraded to the XLP-3000 Syringe)</p> <p>3-Way Valve Assy. For the XL-3000 Syringe Units $<$ 020701SDS</p>	<p>SP6376</p> <p>SP6115</p>
	<p>3-Way Valve Assembly for the Loop (can be used with either the XL-3000 or the XLP-3000; requires minor assembly)</p>	<p>SP6374</p>
	<p>Diluent Bottle – 2 liter</p>	<p>SP6435</p>
	<p>Reagent Bottle Kit – 125mL</p>	<p>SP6436</p>