



**TELEDYNE**  
**CETAC TECHNOLOGIES**  
Everywhereyoulook™



14306 Industrial Road, Omaha, NE 68144, USA  
Tel: 1.800.369.2822 or 1.402.733.2829; Fax: 1.402.733.5292  
Email: sales@cetac.com; Web: www.cetac.com

## **DS-5 MicroFlow Nebulizer**

### **ACCESSORIES & SUPPLIES**

#### **STANDARD VERSION**

**February 4, 2015**

<b>Part Number</b>	<b>Description</b>
SP6030	DS-5 Microflow Nebulizer
SP8157	DS-5 Spray Chamber, Straight Geometry with 12/5 Glass Socket Adapter
SP6024	DS-5 Nebulizer o-ring kit (12 each)
SP5772K	Gas Line Kit, for Agilent 7500/7700/7900/8800 Series ICP-MS, Thermo X-Series 2 and iCAP Q ICP-MS, and Nu Instruments ICP-MS
SP5772L	Gas Line Kit, for PerkinElmer ELAN Series and NexION ICP-MS
SP5772M	Gas Line Kit, for Thermo Element2 & Neptune ICP-MS and Analytik Jena (ex. Bruker, Varian) (All Models) ICP-MS
SP8188	DS-5 Spray Chamber, CEI-100 Style with 12/5 Glass Socket Adapter
SP8241	Torch Adapter for Thermo iCAP Q ICP-MS