



14306 Industrial Road, Omaha, NE 68144, USA
 Tel: 1.800.369.2822 or 1.402.733.2829; Fax: 1.402.733.5292
 Email: cetacsales@teledyne.com; Web: www.cetac.com

QuickTrace™ M-7500

ACCESSORIES & SUPPLIES

STANDARD VERSION

April 29, 2015

Part Number	Description
SP5341	CABLE, RS232 SERIAL
SP5558	SWABS, FOAM TIP
SP5577	THERMAL SAFETY SWITCH; M6000, M7500
SP5580	FAN FILTER ASSEMBLY; M6000, M7500, M8000
SP5582	GLASS CELLS, CELLS ONLY; M6000, M6100, M7500, M7600
SP5584	FLOW RESTRICTORS, (SET OF 2); M6000, M7500
SP5585	OVEN ASSEMBLY; M6000, M7500
SP5586	GLASS CELLS, COMPLETE; M6000, M6100, M7500, M7600
SP5593	BOTTLE KIT, WASTE; M6000, M6100, M7500, M7600, M8000
SP5594	BOTTLE KIT, STANNOUS CHLORIDE REAGENT; M6000, M6100, M7500, M7600,
SP5595	BOTTLE KIT, RINSE; M6000, M6100, M7500, M7600, M8000
SP5603	MERCURY LAMP; M6000, M6100, M7500, M7600
SP5606	CELL O-RING KIT; M6000, M6100, M7500, M7600
SP5608	CELL WINDOW KIT; M6000, M6100, M7500, M7600
SP5612	GLASS BALLAST TANK; M6000, M7500
SP5613	MERCURY LAMP POWER MANAGER; M7500, M7600
SP5614	MERCURY LAMP HOUSING ASSEMBLY; M6000, M6100, M7500, M7600
SP5615	EOFM BOARD; M6000, M6100, M7500, M7600
SP5653	TRAP, MERCURY; M6000, M6100, M7500, M7600, M8000
SP5703	OVEN RADIATOR / HEATER; M6000, M7500
SP5705A	PUMP TUBING, SAMPLE / DRAIN TUBING, YEL-YEL, (PACK OF 12)
SP5705B	PUMP TUBING, STANNOUS CHLORIDE REAGENT TUBING, BLK-BLK, (PACK OF 12)
SP5706	WASTE TUBING, TWO EACH; M6000, M6100, M7500, M7600, M8000
SP5710	GAS FILTER, EXTERNAL, TWO MICRON; M6000, M6100, M7500, M7600, M8000
SP5774	GLS / GLASS BALLAST MOUNT, (2 EACH); M6000, M6100, (M7500 BALLAST MOUNT)
SP5812	GLS CAPILLARIES, (PACK OF 3); M6000, M6100, M7500, M8000
SP5894	NAFION CARTRIDGE REPLACEMENT; M6000, M6100, M7500, M7600, M8000
SP5949	CELL GAS TUBING KIT; M6000, M7500
SP5954	CAMERA 254 Hg FILTERS (SET OF 2); M6000, M6100, M7500, M7600
SP5955	OVEN CONTROLLER; M6000, M7500
SP6033	SAMPLE PROBE, 1.0MM I.D.; M6000, M6100, M7500, M8000
SP6034	CONSUMABLES KIT, 1YR; M7500
SP6035	CELL GAS FITTINGS; M6000, M6100, M7500, M7600

QuickTrace M-7500 Accessories Supplies

SP6036	PERISTALTIC PUMP FITTINGS; M6000, M6100, M7500, M7600
SP6037	GAS EXHAUST FITTINGS (PACK OF 2); M6000, M6100, M7500, M7600
SP6041	GAS LIQUID SEPARATOR; M7500, M8000
SP6042	GLS DRAIN TUBE (PACK OF 3); M6000, M6100, M7500, M7600, M8000
SP6047	FOCAL LENS BLOCK ASSEMBLY(EOFM NOT INCLUDED); M6000, M7500
SP6048	CELL FOCAL LENS, 75mm; M6000, M6100, M7500, M7600
SP6049	MERCURY TRAP MOUNTS, SPRING CLIPS KIT; M6000, M6100, M7500, M7600,
SP6050	DRYER MOUNTS, SPRING CLIPS KIT; M6000 M6100, M7500, M7600, M8000
SP6051	CAMERA 50mm FOCAL LENS (SET OF 2); M6000, M6100, M7500, M7600
SP6052	CAMERA BLOCK; M6000, M6100, M7500, M7600
SP6061	ACQUISITION BOARD; M7500
SP6062	PERISTALTIC PUMP CONTROL BOARD; M6100, M7500
SP6063	INTERCONNECT BOARD; M7500
SP6065	POWER SUPPLY BOARD; M7500
SP6066	VOLTAGE SELECT BOARD; M7500
SP6067	FUSE KIT, 220 V; M7500
SP6068	FUSE KIT, 110V; M7500
SP6069	PRESSURE REGULATOR, INTERNAL; M7500
SP6070	PERISTALTIC PUMP; M7500
SP6071	FAN ASSEMBLY; M7500
SP6072	CAMERA ASSEMBLY; M7500
SP6073	GLS MOUNT / OVERFLOW SENSOR BOARD; M7500, M8000
SP6074	PRESSURE RELIEF VALVE, INTERNAL; M7500
SP6075	POWER SUPPLY, INTERNAL 24V; M7500
SP6076	FIRMWARE, CPU BOARD, VER 2.53 QTA; M7500
SP6077	PROCESSOR BOARD; M7500
SP6079	DETECTOR; M7500
SP6096	COMPLETION KIT, METHOD 1631 / 245.7; M7500, M7600
SP6097	PUMP TUBING, 1631 / 245.7 HYDROXYLAMINE TUBING ORANGE-BLU, 12 EACH
SP6098	CAPILLARY KIT, 1631 / 245.7 HYDROXYLAMINE; M7500, M7600, M8000
SP6155	NAFION ECONOMY PACK - 5 CARTRIDGES; M6000, M6100, M7500, M7600, M8000
SP7013	PERISTALTIC PUMP, STEPPER MOTOR ASSEMBLY; M6100, M7500, M7600, M8000
SP7070	MIXING TEE, METHOD 1631 / 245.7; M7500, M7600
SP7071	MIXING TEE FITTINGS, PEEK, METHOD 1631 / 245.7; M7500, M7600, M8000
SP7072	TUBING, METHOD 1631 / 245.7 HYDROXYLAMINE / SAMPLE TEE TUBE KIT; M7500, M7600
SP7073	TUBING, METHOD 1631 / 245.7 SAMPLE TEE TUBE KIT; M7500, M7600
SP8033	CABLE, PUMP INTERFACE; M7500
SP8034	CABLE, ASX INTERFACE; M7500, M8000
SP8043	PERISTALTIC PUMP MOTOR DRIVE COUPLER; M6100, M7500, M7600, M8000
SP8045	PUMP TUBING, 1631 / 245.7 DRAIN PURPLE / WHITE, 12 EACH; M7500, M7600
SP8049	CAPILLARY, STANNOUS CHLORIDE; M6000, M6100, M7500
SP8062	CAPILLARY REAGENT (10 each), 1631 / 245.7 HYDROXYLAMINE; M7500, M7600
SP8067	CAPILLARY MIXING TEE (10 each), 1631 / 245.7 HYDROXYLAMINE; M7500, M7600
SP8072	CPU / ACQUISITION INTERFACE BOARD; M7500
SP8096	DATA ACQUISITION CABLE; M7500
SP8102	PERISTALTIC PUMP POWER SUPPLY; M6100, M7500
SP8106	PERISTALTIC PUMP MOTOR MOUNT KIT; M6100, M7500, M7600, M8000
SP8108	ENC 500 TRANSFER TUBING

QuickTrace M-7500 Accessories Supplies

SP8153	PUMP TUBING, HGA-100, ORANGE-WHITE, 12 EACH
SP8160	PERISTALTIC PUMP ROLLER HEAD; M6100, M7500, M7600, M8000
SP8161	PERISTALTIC PUMP PRESSURE BRACKET; M6100, M7500, M7600, M8000
SP8162	PERISTALTIC PUMP PRESSURE LEVER; M6100, M7500, M7600, M8000
SP8163	GAS BULKHEAD FITTING; M7500 EXHAUST, M8000 INLET
SP8168	PERISTALTIC PUMP PRESSURE BRACKET SHAFT KIT; M6100, M7500, M7600, M8000
SP8169	PERISTALTIC PUMP PRESSURE ADJUSTMENT LEVER SHAFT KIT; M6100, M7500
SP8187	PRESSURE REGULATOR COMMUNICATION CABLE; M7500
SP8220	PERISTALTIC PUMP, DRIVER BOARD; M7500
SP8252	PERISTALTIC PUMP MOTOR DRIVE BEARING W/ ROLLER SHAFT; M6100, M7500
SP8271	CD S/W QUICKTRACE VER 2.0.4.