


MVX-7100 μ L Workstation is now compatible with Agilent ICP-MS MassHunter Software



ICP-MS MassHunter, the software powering all current Agilent ICP-MS instruments, can now benefit from a plug-in to directly control the Teledyne CETAC MVX-7100 μ L Workstation, bringing micro-sampling capabilities to the Agilent ICP-MS and ICP-QQQ systems!

- Seamlessly operate the entire system via a single dashboard
- Expand ICP-MS capabilities to include:

Micro Sampling Septum Piercing Temperature Control
96 and 384 Well Plate Sampling

- Compatible with the latest MassHunter 4.3 (C.01.03) with Patch 3