

CETAC TECHNOLOGIES

SERVICE BULLETIN

| | |
|--|------------------------------|
| Subject: Automation systems pumps | Bulletin No: SB1012 |
| | Date Issued: 12/31/07 |
| | Priority: Low |

Product: ADX, ASX-110FR, ASX-112FR, ASX-130, ASX-260, ASX-410, ASX-520, ASX-520HS, ASX-750, EXR-8, SDS-550

Serial Numbers: >080101XXXXXXXX
Serial Numbers start by month (MM), year (YY) and units built that month (i.e., MMYyxx)

Service Location: Customer Site CETAC Technologies

Responsibility: Customer Trained CETAC Representative

Replacement Parts: Customer CETAC Technologies

DESCRIPTION

Due to design changes to the products listed above CETAC Technologies has a new pump for the Automation line of products.

IMPLEMENTATION

Replace the peri pump with the new SP7002A following the direction in the documentation provided.

The pumps are replaced on failure and are a direct replacement for SP7002. This is a consumable part.

The new pumps will be available after January 3, 2008

PARTS

| Part Number | Description |
|-------------|--------------------|
| SP7002A | Pump Motor Kit, WM |

AVAILABLE DOCUMENTATION