

Dust Cover Enclosures for Autosamplers



ACCESSORIES AND PARTS CATALOG FOR THE ENC-560DC, ENC-280DC, AND ENC-860DC DUST COVER ENCLOSURES

Dust Cover Accessories and Parts



Enclosure Fan—positive air flow for use with HEPA filter when vent is not available. Provides a flow rate of 121 CFM.

North America	32-0099-045
Europe	32-0100-045
Australia	32-0101-045



HEPA Filter (Requires adapter plate)

Filter with adapter plate 32-0162-045



Filter replacement only 32-0102-045



Exhaust Flange,

2.5" Diameter	32-0288-045
3.0" Diameter	32-0289-045
3.5" Diameter	32-0290-045
4.0" Diameter	32-0103-045



Autosampler Base

For ASX-280 (ENC-280DC)	32-0105-045
For ASX-560 (ENC-560DC)	32-0104-045
For XLR-860 (ENC-860DC)	32-0106-045

Ordering

To order, contact:

Phone (Toll-Free, USA): 800-369-2822
Phone (Outside USA): +1-402-733-2829
Fax: +1-402-733-5292
e-mail: cetacsales@teledyne.com
www.cetac.com

Teledyne CETAC Technologies
14306 Industrial Road
Omaha, Nebraska 68144, USA

For questions about service and repair, call us at the numbers above or e-mail
cetacservice@teledyne.com