

# APS-7650G GRAVIMETRIC SAMPLE PREP STATION



## ACCESSORIES AND PARTS CATALOG

Catalog part number 32-1265-062 Rev 00, January 2022

# Accessories and Supplies



**Probe Assembly**

32-1266-047



**Stainless Steel Tube**

32-1057-030

18" Tube



**Gripper Tips**

32-1242-044



**25 mL Syringe**

32-1054-031



**2.5 mL Syringe with Pipettor**

32-1203-047



**3-Way Valve**

SP6376

## APS-7650G ACCESSORIES AND PARTS



Polyurethane Tubing

SP6300




Waste Station

32-1196-047



Barbed Fitting Kit

32-1244-044




Rinse/Drain Tubing Hookup Kit

Material	I.D.	Quantity	Part Number
Tygon® Fuel and Lubricant	1/8"	2 x 84"	SP6253
Superthane®	1/8"	2 x 84"	SP6431
Tygon® Fuel and Lubricant	3/16"	2 x 72"	SP6254
Superthane®	3/16"	2 x 72"	SP6432



45 Position Bottle Rack

SP6274



90 Cup Rack

SP6276

## APS-7650G ACCESSORIES AND PARTS



**60 Position Metal Tube Rack**  
Holds 15 mL tubes

SP7163



**80 Position Metal Tube Rack**  
Holds 15 mL tubes

SP7007



**Sample Cup, 2 mL**  
500 pieces

SP5178P



**15 mL (17mm x 100mm) sample vial**  
1000 pieces

SP5178M



**4 oz Bottle with Cap**

SP5178S



**Pipette Tips – 1 Box (1000 tips)**

32-1221-044

## APS-7650G ACCESSORIES AND PARTS



**Pipette Waste Bin**

32-1127-045

### Power Cords



Where Used	Description	Part Number
USA, Canada, Japan, Mexico	Type B NEMA 5-15 IEC 60906-2 100-127 V	SP5003
Europe, Russia	Type F CEE 7/7 220-240 V	SP5004



**USB Cable**

SP6375



**RS-232 Interface Cable**

3 Meters

SP7441



**6-Pin DIN Cable**

32-1106-036

# Repair Parts




**Gripper Module**

32-1123-047



**Syringe Motor and Level Sensor Assembly**

32-1202-047



**APS-7650G Z-Drive Assembly**

32-1243-047



**APS-7650G Vortex Module**

32-1194-047



**Precision Balance**

32-1195-047

# APS-7650G ACCESSORIES AND PARTS

**APS-7650G I/O Board**

32-1201-047



**APS-7650G Main Board**

32-1245-047



**Sidekick Board**

SP7138



**W-Axis Board**

Includes two mounting screws

SP7169



**APS-7650G Power Supply**

32-1198-047



## APS-7650G ACCESSORIES AND PARTS



**X-Axis Drive Belt**

SP5222



**Balance Cover**

32-1093-037



**SPM-700G**  
(Syringe not included)

32-1169-047



## APS-7650G ACCESSORIES AND PARTS

# Ordering

To order, contact:

Phone (Toll-Free, USA): 800-369-2822  
Phone (Outside USA): +1-402-733-2829  
e-mail: [cetacsales@teledyne.com](mailto:cetacsales@teledyne.com)  
[www.teledynecetac.com](http://www.teledynecetac.com)

Teledyne CETAC Technologies  
14306 Industrial Road  
Omaha, Nebraska 68144, USA

For questions about service and repair, call us at the numbers above or e-mail  
[cetacservice@teledyne.com](mailto:cetacservice@teledyne.com)

Trademark Acknowledgements: PharMed and Tygon are registered trademarks of Saint-Gobain Performance Plastics. Viton® is a registered trademark of E.I. du Pont de Nemours and Company. Superthane is a registered trademark of NewAge Industries, Inc. KYNAR® is a registered trademark of Arkema Inc.