

U5000AT+ & U6000AT+ Essential Spare Parts

Avoid downtime with an on-hand supply of these recommended parts.

One-Year Consumable Kit

part number
SP7679

Includes:

- 1 Transducer Assembly
- 2 each 2.5A, 5A, and 3A Fuses
- 1 set of Drain Pump Tubing
- 1 Argon In flow restrictor
- 2 ea. PEEK sample inlet tubing
- 60 in. (1.5m) PVC Sample out tubing, 3/16" ID
- 180 in. (4.6m) PVC waste tubing, 1/8" ID
- 12 pk peristaltic pump tubing, 1.3 mm ID
- 24 in. (0.6m) of drain pump tubing
- 1 spare aerosol chamber o-ring
- 36 in. (0.9m) 0.8 mm ID PTFE sample tubing
- 24 in. (0.6m) Rinse tubing 3/32" ID

**1 Year
Kit**



**Order a
spare**



**Fuse Kit
SP5410**

**6
months**



**Transducer
O-ring Kit**

Two o-ring
set. 6 month
recommended
replacement

SP5168B

**6-12
months**



**Spare Part
Kit**

6-12 month
recommended
replacement

SP5447

**6
months**



**Transducer
Assembly**

Includes o-ring.
6 month
recommended
replacement

SP5125

**6
months**




**PEEK
Sample Inlet
Tubing**

6 month
recommended
replacement

SP5150B

**6
months**



**Drain Pump
Tubing and
Connector
Kit**

6 month
recommended
replacement

SP5366

Order Now: ☎ 402.733.2829 ✉ cetacsales@teledyne.com 🖨 store.teledynecetac.com

Useful Accessories for the U5000AT+ & U6000AT+

Add convenience to your laboratory workflow with these helpful items.



Utility Cart

Provides a mobile platform for the U5000AT+ or U6000AT+ ultrasonic nebulizers and spare parts.

Part Number

450492



Waste Bottle Kit

Large capacity (15 Liter) LDPE bottle for liquid sample waste.

Part Number

SP8320




Extension Tubing Kit

Allows easy extension of the U5000AT+ or U6000AT+ sample inlet tubing to a peristaltic pump.

Part Number

SP5154



Cutting Block With Blade

Guarantees an accurate, angled cut for the U5000AT+ or U6000AT+ sample inlet tubing.

Part Number

SP5170

Order Now: ☎ 402.733.2829 ✉ cetacsales@teledyne.com 🌐 store.teledynecetac.com

Teledyne CETAC Technologies, 14306 Industrial Road, Omaha, NE 68144 USA
+1 402.733.2829 | +1 800.369.2822 | cetacsales@teledyne.com | www.teledynecetac.com