

Useful Accessories for the Teledyne CETAC Ultrasonic Nebulizers U5000AT+ or U6000AT+



Utility cart

Provides a mobile platform for the U5000AT+ or U6000AT+ ultrasonic nebulizers and spare parts.


SP5158



Waste bottle kit

Large capacity (10 Liter) bottle for liquid sample waste. Bottle lid is equipped with three ports for the U5000AT+ or U6000AT+ drain lines and one vent port.

SP7039



Cutting block with blade

Guarantees an accurate, angled cut for the U5000AT+ or U6000AT+ sample inlet tubing.

SP5170



Extension tubing kit

Allows easy extension of the U5000AT+ or U6000AT+ sample inlet tubing to a peristaltic pump.

SP5154

Please contact Teledyne CETAC (CETACsales@Teledyne.com) for ordering information.