



ASX-110FR AND ASX-112FR MICRO AUTOSAMPLER ACCESSORIES AND SUPPLIES CATALOG



Teledyne CETAC Technologies
14306 Industrial Road
Omaha, Nebraska 68144, USA
Phone: (402) 733-2829
Fax: (402) 733-5292
e-mail: custserv@cetac.com
www.cetac.com
620083 Rev1, May, 2014

ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	<p>Sample Probe</p> <p>0.010" (0.3mm) ID Carbon – Black Band</p> <p>0.020" (0.5mm) ID – Blue Band</p> <p>0.035" (0.9mm) ID Carbon – Yellow Band</p>	<p>SP6328</p> <p>SP6397</p> <p>SP6329</p>
	<p>Short Rack Kit</p> <p>24 Position Rack with Cover</p> <p>With qty 24 - 2.0mL Polypropylene Vials</p> <p>With qty 24 – 1.5mL PFA Vials</p>	<p>SP6321</p> <p>SP6322</p>
	<p>Short Rack Kit</p> <p>48 Position Rack</p> <p>Qty 48 - 0.5mL Polypropylene Vials, with Screw On Caps</p>	<p>SP6323</p>
	<p>Short Rack Kit</p> <p>48 Position Rack with Cover</p> <p>Qty 48 - 0.5mL Teflon Vials</p>	<p>SP6324</p>

ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	<p>Short Rack Kit 96 Position Rack with Cover Qty 96 - 1.0mL Polypropylene Vials</p>	<p>SP6325</p>
	<p>Short Rack Tray Adaptor Kit</p>	<p>SP6337</p>
	<p>Tall Rack Kit 30 Position, Bel-Art Rack Qty 30 - 14mL Polypropylene Vials</p>	<p>SP6326</p>
	<p>Tall Rack Kit 42 Position, Bel-Art Rack Qty 42 - 7mL Polypropylene Vials</p>	<p>SP6327</p>



ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	Rotary Tray	SP6336
	Sample Vials (for use with 24 position short racks) 2.0mL Polypropylene Vials (qty 1,000) 2.0mL Polypropylene Vials (qty 100) 1.5mL PFA Vials (qty 24)	SP5540 SP5540A SP5472
	Sample Vials (for use with 48 position short racks) 0.5mL Polypropylene Vials (qty 1,000) Screw Cap for 0.5mL Polypropylene Vial (qty 1,000)	SP5479 SP5511
	Sample Vials (for use with 48 position short racks) 0.5mL Teflon Vials (qty 48)	SP5477


ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	<p>Sample Vials (for use with 96 position short racks)</p> <p>1.0mL Polypropylene vials (qty 1,000)</p>	<p>SP5480</p>
	<p>Sample Vials (use with Bel-Art racks only)</p> <p>14mL Polypropylene, 16mm x 100mm (1000/kit) SP5178B 14mL Polypropylene, 16mm x 100mm (5000/kit) SP5178B5 14mL Polypropylene, 16mm x 100mm (250/kit) SP5178F</p> <p>8mL Polypropylene, 13mm x 100mm (1000/kit) SP5178K 8mL Polypropylene, 13mm x 100mm (250/kit) SP5178L</p> <p>(Round bottom without cap)</p> <p>7mL Polypropylene, 13mm x 82mm (1000/kit) SP5178A 7mL Polypropylene, 13mm x 82mm (250/kit) SP5178E</p> <p>(Flat bottom without cap)</p>	
	<p>Short Standards Vials (use with short racks only)</p> <p>4mL PFA Standards Vials (qty 9)</p>	<p>SP6332</p>
	<p>Short Standards Vials (use with short racks only)</p> <p>20mL PFA Standards Vials (qty 5)</p>	<p>SP6333</p>
	<p>Tall Standards Vials (use with tall racks only)</p> <p>10mL PFA Standards Vials (qty 9)</p>	<p>SP6334</p>

ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	<p>Tall Standards Vials (use with tall racks only)</p> <p>30mL PFA Standards Vials (qty 5)</p>	<p>SP6335</p>
	<p>Short Rinse Kit for Peri Pump</p> <p>With Polypropylene Rinse Station (ASX-110FR)</p> <p>With Polypropylene Rinse Station (ASX-112FR)</p>	<p>SP6315</p> <p>SP6353</p>
	<p>Tall Rinse Kit for Peri Pump</p> <p>With Polypropylene Rinse Station (ASX-110FR)</p> <p>With Polypropylene Rinse Station (ASX-112FR)</p>	<p>SP6316</p> <p>SP6354</p>
	<p>Short Rinse Kit for Gas Displacement Models</p> <p>With Polypropylene Rinse Station (ASX-110FR)</p> <p>With PFA Rinse Station (ASX-110FR)</p> <p>With Polypropylene Rinse Station (ASX-112FR)</p> <p>With PFA Rinse Station (ASX-112FR)</p>	<p>SP6317</p> <p>SP6319</p> <p>SP6355</p> <p>SP6357</p>
	<p>Tall Rinse Kit for Gas Displacement Models</p> <p>With Polypropylene Rinse Station (ASX-110FR)</p> <p>With PFA Rinse Station (ASX-110FR)</p> <p>With Polypropylene Rinse Station (ASX-112FR)</p> <p>With PFA Rinse Station (ASX-112FR)</p>	<p>SP6318</p> <p>SP6320</p> <p>SP6356</p> <p>SP6358</p>

ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	Peri Pump Motor Kit 2-Channel	SP7002A
	Switch / Power Input (6 pin input)	SP6104
	Main Ribbon Cable	SP6690
	Peri Pump Rinse Tubing Kit (ASX-110FR) (Tygon) Peri Pump Rinse Tubing Kit (ASX-112FR) (Tygon)	SP6330 SP6359

ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	<p>Gas Displacement Rinse Tubing Kit (ASX-110FR) (Tygon / PFA / PTFE)</p> <p>Gas Displacement Rinse Tubing Kit (ASX-112FR) (Tygon / PFA / PTFE)</p>	<p>SP6331</p> <p>SP6360</p>
	<p>ASX-110FR Cover</p>	<p>SP6361</p>
	<p>ASX-112FR Cover</p>	<p>SP6362</p>
	<p>ASX-110 / ASX-112 Vent Port Kit</p>	<p>SP6876</p>





ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	Bottle, Gas Displacement Reservoir, PFA	SP6341
	Z-Drive PFA Tubing Replacement Kit	SP6338
	Fastener Service Kit	SP6621

ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	Rinse Station Thumb Screws (qty 2)	SP6363
	Probe Holder Thumb Screw Kit (Qty 4) (Qty 6)	SP6301 SP7089
	Y-Arm Attachment Screw (qty 2)	SP6364
	24 Volt Power Supply	SP6103
	Serial Cable	SP5341

ASX-110FR & ASX-112FR AUTO SAMPLER SPARE PARTS

IMAGE	DESCRIPTION	PART #
	ASX-110FR, ASX-112FR Main Board	SP6339
	ASX-110FR, ASX-112FR Rabbit Board with Firmware	SP6340
	Drain Pump Tubing and Connector Kit (PharMed L/S 14)	SP6089
	Drain Pump Tubing and Connector Kit (PharMed L/S 13 – Low Flowrate Option)	SP6190

PharMed is a registered trademark of Norton Co. Tygon is a registered trademark of Norton Co