





DSA-7
DISCRETE SAMPLING ACCESSORY
ACCESSORIES AND SUPPLIES CATALOG



CETAC Technologies
14306 Industrial Road
Omaha, Nebraska 68144, USA
Phone: (402) 733-2829
Fax: (402) 733-5292
e-mail: cetacsales@teledyne.com
www.teledynecetac.com
620084 August, 2018

Tygon and Pharmed is a registered trademark of.
Saint-Gobain Performance Plastics.
Tefzel is a registered trademark of DuPont

DSA-7 Discrete Sampling Accessory Spare Parts Catalog

IMAGE	DESCRIPTION	PART #
	<p>Sample Probe, Carbon Fiber, 1.0mm (Standard with DSA)</p>	<p>SP6408C</p>
	<p>Sample Probe, Ultem, 1.0mm (Optional with DSA-7)</p>	<p>SP6405</p>
	<p>Sample Probe, Carbon Fiber, 0.5mm (Standard with DSA. Used with Autosampler when not utilizing DSA-7)</p>	<p>SP5800C</p>
	<p>Sample Loop Assy.—0.70mL Sample Loop Assy.—1.00mL Sample Loop Assy.—1.25mL Sample Loop Assy.—1.50mL Sample Loop Assy.—2.00mL Sample Loop Assy.—2.50mL Sample Loop Assy.—3.00mL Sample Loop Assy.—3.50mL Sample Loop Assy.—4.00mL Sample Loop Assy.—4.50mL Sample Loop Assy.—5.00mL Sample Loop Assy.—5.25mL</p>	<p>SP6415 SP6384 SP6385 SP6386 SP6387 SP6388 SP6389 SP6390 SP6391 SP6392 SP6393 SP6394</p>

DSA-7 Discrete Sampling Accessory Spare Parts Catalog

IMAGE	DESCRIPTION	PART #
	<p>Barbed Fitting Kit, ¼-28 male to 1/8" barb (qty 6) (Used on serial number < 100901DSA7)</p>	<p>SP6298</p>
	<p>Nut/Ferrule & Fitting Kit (Two green flangeless connectors with 1/16" ferrules) (Two blue flangeless connectors with 1/16" ferrules) (Four Tefzel flangeless connectors with 1/16" ferrules) (Two ¼"-28 to 1/8" barb fittings)</p>	<p>SP6399</p>
	<p>6-Port Valve PPS (polyphenylene sulphide) construction (Used on serial number < 100901DSA7)</p>	<p>SP6400</p>
	<p>6-Port Valve (S.B.) PPS (polyphenylene sulphide) construction (Used on serial number ≥ 100901DSA7)</p>	<p>SP6541</p>





DSA-7 Discrete Sampling Accessory Spare Parts Catalog

IMAGE	DESCRIPTION	PART #
	<p>6-Port Valve Coupler</p>	<p>SP6412</p>
	<p>Rinse Tubing Kit (Tygon/Teflon/PharMed)</p>	<p>SP6406</p>
	<p>Rinse Tubing Kit, Ultra-Clean (Used for ultra-clean applications) (Teflon lined Tygon/Teflon/PharMed)</p>	<p>SP6407</p>




DSA-7 Discrete Sampling Accessory Spare Parts Catalog

IMAGE	DESCRIPTION	PART #
	Carrier/Rinse Bottle	SP6409
	Vacuum Pump	SP6404
	3-Port Anti-Siphon Valve (Used on serial number < 030901DSA7)	SP6403
	2-Port Anti-Siphon Valve (Used on s/n \geq 030901DSA7 < 100901DSA7)	SP6473

DSA-7 Discrete Sampling Accessory Spare Parts Catalog

IMAGE	DESCRIPTION	PART #
	<p>Daughter Board Kit (Installs in ASX-520/520HS/EXR-8 to allow interface with DSA-7)</p>	<p>SP6147</p>
	<p>Daughter Board Kit (Installs in ASX-260 to allow interface with DSA-7)</p>	<p>SP6127</p>
	<p>DSA-7 Enabled Rabbit Processor Module (CPU) with Firmware (for Standard Speed Autosamplers)</p> <p>Used on ASX-260, ASX520, EXR-8 <u>Standard Speed Autosamplers</u> with DSA-7.</p> <p><i>(Used on standard speed autosamplers connecting to DSA-7 with serial numbers \geq 100901DSA7)</i></p>	<p>SP6543</p>
	<p>DSA-7 Enabled Rabbit Processor Module (CPU) with Firmware (for High Speed Autosamplers)</p> <p>Used on ASX520HS and EXR-8 <u>High Speed Autosamplers</u> with DSA-7.</p> <p><i>(Used on high speed autosamplers connecting to DSA-7 with serial numbers \geq 100901DSA7)</i></p>	<p>SP6544</p>

DSA-7 Discrete Sampling Accessory Spare Parts Catalog

IMAGE	DESCRIPTION	PART #
	<p>Vacuum Pump Flow Restrictor Kit (Used w/2-port Anti-Siphon Valve) (Used on s/n \geq 030901DSA7 < 100901DSA7)</p>	<p>SP6474</p>
	<p>Vacuum Pump Flow Restrictor Line (S.B.) (Used on s/n \geq 100901DSA7)</p>	<p>SP6527</p>
	<p>Carrier Solution Tee Assembly (S.B.) (Used on s/n \geq 100901DSA7)</p>	<p>SP6528</p>