


Pump Up Protein Purification

An off-the-shelf solution for automated high throughput sample injection system with the ÄKTA™ pure chromatography system

The ASX-560 Autosampler provides an off-the-shelf solution for increasing sample throughput while decreasing hands-on time with your preparative LC instrument, increasing the amount of time available for your research.

The ASX-560 Autosampler allows for process consistency as well as unattended purification of:

- Proteins
- Enzymes
- Other biomolecules



- Off-the-shelf solution for fully automated protein purification
- Manage sample volumes from 2 mL to up to 250 mL
- Highly reproducible recoveries
- No modifications to the ÄKTA™ system
- Increase throughput while decreasing hands-on time




TELEDYNE
CETAC TECHNOLOGIES
Everywhereyoulook™

ASX-560 Autosampler Options

Flexibility

Manage sample volumes of 8 mL up to 250 mL. Multiple rack options are available, including adapters for large-volume bottles. The ASX-560 Autosampler is also compatible with other standard, ready-made sample racks.



25 Position 250 mL Rack Adapter SP7640



18/21 Position Adjustable Rack 32-0748-047

Protection

An optional clean enclosure guards against environmental impurities to protect both operator and equipment.



ENC-560 DC enclosure shown with ASX-560 (not included)

ENC-560 DC for ASX-560 ENC-99-0010

HEPA filter with adapter plate 32-0162-045

Ordering

ASX-560 for ÄKTA pure	A56-99-0026
ASX-280 for ÄKTA pure	A28-99-0026

Software

Universal AScript autosampler scripting software syncs with the ÄKTA system via COM cable.

Communication Cables

Cable for ÄKTA pure	SP8318
---------------------	--------

Capacities

Volume	ASX-560	
	Racks	Samples
10 mL	4 x 90	360
15 mL	4 x 60	240
50 mL	4 x 21	84
250 mL	1 x 18	18
250 mL	1 x 21	21
250 mL	1 x 25	25