

**Certificate of Decontamination**

This form is for use when returning equipment to CETAC Technologies for repair, refurbishing, replacement or from demonstration. Any equipment that has been located in a biological or chemical laboratory must be assessed and, if necessary, decontaminated prior to acceptance by CETAC. CETAC Technologies **does not accept** equipment that has been exposed to **radioactive contaminants** or **BioSafety level 4 laboratories**. Please complete this form, sign it and send a copy of the form with the instrument, placed on the outside of the package.

**Company** \_\_\_\_\_

**RMA#** \_\_\_\_\_

**Instrument** \_\_\_\_\_

**Reason for return** \_\_\_\_\_

**Decontamination information:**

Has the instrument been exposed to any infectious agents/biological samples assigned to biosafety levels 2 or 3? \_\_\_\_\_

If so, please indicate the agent and bio-safety hazard class. \_\_\_\_\_  
(Please note CETAC Technologies does not accept equipment removed from biosafety level 4 labs)

Has the instrument been exposed to chemical solutions? \_\_\_\_\_

If so, please list the chemicals used in the instrument with their concentration levels. \_\_\_\_\_

**Decontamination Procedure:**

Please indicate how the instrument was disinfected or cleaned by checking one or more of the decontamination procedures described below:

**Equipment not contaminated** Customer/user certifies that the returned equipment has not come into contact with, or contain any infectious agents or chemical solutions.

**Contact with infectious agents** The instrument was thoroughly decontaminated by spraying it or wiping down with a 1:10 dilution of Clorox brand bleach, or 1/2% solution of Lysol IC brand disinfectant, or equivalent. The disinfectant remained on the unit for a minimum of 10 minutes before flushing with water.

**Chemical contamination** The equipment has been thoroughly decontaminated by wiping and rinsing areas associated with the chemicals with an ammonia based cleaner such as Windex with ammonia and then flushed with di water.

By accepting authorization to return the product, I assume all responsibility for all chemical and biological decontamination and cleaning. I understand that CETAC Technologies has no obligation to repair or service any product unless adequate information is provided to ensure the safety of all personnel handling the exposed goods.

Printed Name: \_\_\_\_\_

Company Name: \_\_\_\_\_

Phone Number: \_\_\_\_\_

Signature: \_\_\_\_\_ Date: \_\_\_\_\_